FLOWOX™

Enhancement of blood flow to lower limbs

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


INSTRUCTION MANUAL [EN]

WELCOME

FlowOx™ is a medical device in accordance with the Medical Device Regulation 2017/745, which improves blood flow to your lower limb. The patented technology is based on several years of research and development.

Please carefully read this instruction manual before using the device. The Quick Guide [0] should be kept handy.

Otivio and its distributors decline any responsibility for any use of FlowOx™ that does not comply with this instruction manual.

HOW IT WORKS

The user puts the foot and leg into a Pressure Chamber [1] which is connected to the Control Unit [4]. A pump in the Control Unit creates a pre-set cycle of 10 seconds with negative pressure of -40 mmHg +-5 mmHg and 7 seconds with ambient air pressure, both targets are reached within 8 seconds and 3 seconds, respectively. The negative pressure facilitates blood flow down into the limb and to the skin. The Control Unit [4] is a multiple patient multiple use device; hence it may be used for multiple patients during its lifetime. The Pressure Chamber [1] directly contacts the skin; hence it is a single patient multiple use device [59] and should be used by a single patient only.

Use FlowOx™ as recommended by your clinician, normally for 2 hours daily, divided into at least two sessions. Follow up your clinician recommendations for treatment modalities.

INTENDED USE, INDICATIONS FOR USE AND CONTRA-INDICATIONS

Intended Use:

The FlowOx™ System is intended to be used for the treatment of conditions and diseases of the foot, ankle and lower leg that are associated with impaired blood flow. It enhances peripheral blood flow by application of negative pressure pulses to the treated area. It is intended to be used in both a professional and home care environment.

Indications:

FlowOx™ is indicated for the treatment of the following conditions in adults (over 18 years of age):

- Peripheral arterial disease.
- Diabetic and non-diabetic arteriopathic foot and leg ulcers.
- Intermittent claudication and rest pain.
- Foot and lower leg ulcers of mixed aetiology and/or pain/spasticity or discomfort associated with immobility such as paraplegia following spinal cord injury or other neurodegenerative diseases such as multiple sclerosis.

Intended users:

FlowOx™ 2.0 is intended for adult users aged 18 and above. There are several different categories of intended users, each utilizing different aspects of the FlowOx™ system. The primary intended users include:

- Patients: Primarily elderly individuals with limited physical capabilities.
- Clinicians: Medically trained professionals such as doctors or registered nurses.
- Nurse Assistants: Trained caregivers working in nursing homes or providing in-home care.
- Carers: Individuals assisting the patient at home.

These categories encompass the range of users for the FlowOx™ system.

Contraindications and Warnings:

There are no known contraindications for the device.

- Due to an absence of data, caution should be taken in its use in non-adults, under age of 18 years and in pregnant women.
- Due to the close contact of the pressure chamber and the skin, underlying local conditions such as psoriasis, eczema, varicose veins or skin infections may become irritated or exacerbated.
- Reperfusion pain is experienced by some patients.

- Due to increased blood flow to the treated area, ulcers may experience increased exudate production or bleeding. Should these occur, a super absorbent dressing may be required and/or dressings may need to be changed more often or treatment terminated.
- Underlying conditions and co-morbidities should be stabilised and under control prior to treatment to be compatible with the treatment period.
- If any medical conditions become aggravated, treatment should be stopped, and medical advice obtained.
- Some patients may face challenges in tolerating FlowOx treatment, as it necessitates extended periods of sedentary activity, and initial treatment may exacerbate feelings of discomfort, pain, and spasticity. In such instances, it is advisable to discontinue the treatment.

Treatment Period:

FlowOx™ use is recommended for 2 hours per day and clinical improvement in symptoms usually requires a minimum of 8 weeks treatment. Treatment is applied whilst the patient is seated, and it does not interfere with other activities normally enjoyed whilst seated.

INITIAL SETUP OF FLOWOX™ BY CLINICIAN

Please read this FlowOx™ instruction manual and refer to the drawings at the beginning of the instruction manual.

- (A)** Unpack and verify that all items are provided and free from damage (see overview picture). If any parts are missing or damaged, please contact your clinician or distributor. The Control Unit box contains the Quick Guide [0], Instructions for use (IFU) Manual [0], the Control Unit [4], the power supply [2], Timecard [5] (if supplied), USB drive with FlowOx™ Software [32] data transfer cable [11]. The Pressure Chamber box contains the Quick Guide [0], Instructions for use (IFU) Manual [0], the PCH [1], 2x AbsorbGel pouch [26], Positioner [16], Padding [17], 2x Seal [20] and Strap [23]. The Positioner [16], Padding [17], Seal [20], and Strap [23] are applied parts of the system and the materials are non-allergenic.
- (B)** Connect the External Power Supply [2] to the Power Supply Port [3] on the Control Unit [4] and plug into an electrical socket. **Note:** Use of external power supply [2] other than that specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation. The type of External Power Supply [2] is Meanwell GSM60B24-P1L. **Note:** Make sure the External Power Supply [2] can easily be disconnected from the Control Unit [4] during use.
- (C) Check FlowOx™ Control Unit for the timecard:** To check if the Control Unit [4] is activated with a Timecard, press the On/Off Button [7]. A timecard is activated if the Control Unit [4] starts the treatment. If instead, the Control Unit Display [10] shows the text “Oh” and a General Error symbol [35], then no Timecard [5] is yet activated. The FlowOx™ Control Unit [4] can only be operated when it has been activated with a Timecard [5]. To add time, do the following: **|a|** Check that the Control Unit [4] is switched off. **|b|** Insert the Timecard [5] into the Timecard Port [6]. **|c|** Press the On/Off Button [7]. When the Control Unit [4] has read the Timecard [5], an hourglass symbol, along with total time left, will be displayed on the Multi Screen [8] and the Smiley Face symbol [9] will appear on the Display [10]. The Control Unit [4] will shut down automatically. **|d|** Remove the Timecard [5] from the Control Unit [4] and dispose of it immediately. **Note:** Timecards [5] come with varying amounts of hours.
- (D) Set time (hours).** Time is set manually on Control Unit [4] or by synchronizing with a Personal Computer (PC) using the USB Cable [11]. USB cable supplied with the Control Unit is a 2 m long, black type A-B USB cable. **Manual time setup:** **|a|** Push the On/Off Button [7] for 10 seconds. The Clock Set symbol [12] lights up. The hour blinks on the Multi Screen [8] **|b|** Press the On/Off Button [7] repeatedly until the correct local hour is set. **Note:** minutes do not have to be set. **|c|** Wait for 10 seconds. The local time will be briefly displayed. **|d|** The Control Unit [4] shuts down automatically afterwards. **Synchronized PC time setup:** See section “Viewing Treatment Data (Clinician)” for information on how to setup the software **|a|** Connect the USB Cable [11] to the USB Cable Port [13] on the Control Unit [4] and then to the PC **|b|** Start the FlowOx™ Software **|c|** Press the On/Off Button [7] on the Control Unit [4] **|d|** Click on the Cable icon [14] and then the Clock Set icon [15] in the FlowOx™ Software **|e|** The Control Unit [4] adopts the local time from the PC and the Clock Set symbol [12] is lit. The time is shown on the Multi Screen [8] after it is set. The Control Unit [4] shuts down automatically in 15 seconds. **|f|** Disconnect the USB Cable [11] at both ends.

FIRST TIME FITTING

Refer to the drawings from the cover page of this instruction manual. In all the following steps, the clinician should train the patient. **The clinician will assess whether the patient will need someone to help them with using FlowOx™. FlowOx™ is a sophisticated medical device. It may take the patient several training sessions to become familiar enough to use the device routinely.**

- (A) Preparing FlowOx™:** (a) Insert the Positioner [16] into the Pressure Chamber [1] and align its grooves with the ridges on the floor of the Pressure Chamber [1]. The Positioner [16] is adjustable forwards and backwards and for left or right foot usage. The foot shall be placed on the Positioner, so the foot alignment matches the Left Foot symbol [57] or Right Foot symbol [58]. (b) Insert the Padding [17] into the Pressure Chamber [1] so that the Fastening Band [18] is towards the top. Start by pushing its two Air Vents [19] into the holes on the inside of the Pressure Chamber [1]. Attach the four holes on the Fastening Band [18] of the Padding [17] to the four studs on the inside of the Pressure Chamber neck. (c) Place the Seal [20] around the top of the Pressure Chamber [1]. The Seal has a hole in the front section which should be placed on the Stud [21] on the front top of the Pressure Chamber [1]. **Note:** The Seal Seam [22] should be on the outside (d) Place one of the holes of the Strap [23] over the same stud that the Seal [20] is fastened to. Wrap the Strap [23] around the Seal [20] edges and secure it by placing the other Strap [23] hole over the same Stud [21] (e) Connect the Hose with Filter [24] to the Hose Connector [25] on the Control Unit [4]. **Note:** Ensure that the connection is tight. (f) Connect the External Power Supply [2] to the Power Supply Port [3] on the Control Unit [4] and then to an electrical socket. **Note:** No other power supply should be used.
- (B) Fitting FlowOx™.** (a) Check that the Pressure Chamber [1] fits the patient and is comfortable. Only the arch of the foot should rest lightly on the top of the Positioner [16]. **Note:** The foot should not touch the bottom of the Pressure Chamber [1]. The Positioner [16] is not a weight bearing support device. (b) The Seal [20] should be comfortable on the bare calf. (c) Ask the patient to put their Pressure Chamber [1] on and take it off several times. Check that the patient can do this without harming the skin or disrupting any wound dressings. **Note:** It is recommended that a non-shedding sock/tubular bandage is used to cover any wound dressings and the areas of the leg which touch the Padding [17] and Positioner [16]. (d) If the test was OK, the patient can use the device alone as instructed in the section “Daily Use of FlowOx™”. **Note:** training should include: how to set up the device at home, put on and take off the Pressure Chamber [1], operating the system, cleaning, and how to care for the device, adding time with the Timecard [5], and discussion of hazards, e.g. tripping and disposal of small parts.

Please refer to the User Instructions videos on the 'Help' page of the website (www.flowox.com) for further information.

DAILY USE OF FLOWOX™

- **Putting on the Pressure Chamber [1]:** Once the system has been assembled (see “Preparing FlowOx™ above), (a) Roll down the Seal [20] over the Pressure Chamber [1] neck to allow the leg to be put into the opening. (b) While sitting on a chair, insert the foot to be treated into the Pressure Chamber [1]. (c) Make sure that the Positioner [16] is correctly placed for the foot to be treated as advised by the clinician. (d) Roll the Seal [20] up over the leg and make sure it fits snugly around the bare leg. **Note:** Physical contact should only occur between the patient’s leg, the Positioner [16], the Padding [17] and the Seal [20] and caution should be taken to avoid pressure points between the leg and the device.
- **Treatment:** (a) Place the leg into the Pressure Chamber [1] and press the On/Off Button [7] on the Control Unit [4] to start the treatment. (b) Make sure the leg and foot are centrally placed inside the Pressure Chamber [1] and not touching the outside walls during inflation of the Padding [17] (10-15 sec). See image [29]. (c) Treatment begins. Time counts upwards on the Multi Screen [8]. The Treatment Ongoing symbol [27], along with a progress bar (shown first minute only) will appear on the Multi Screen [8]. (d) Place the Pressure Chamber [1] at a comfortable distance from the chair and if needed, change position of the Pressure Chamber [1] during treatment to increase comfort. **Note:** The patient should remain seated during treatment ● After 2 hours a Smiley Face symbol [9] lights up on the Multi Screen [8] ● All treatment use is recorded on the device ● The FlowOx™ continues working for 12 hours unless it is stopped ● If condensation appears on the inside of the Pressure Chamber [1] when in use, insert 1-2 Absorgel Pouches [26] into the Pressure Chamber [1] floor. Note: An Absorgel Pouch [26] should be used for a maximum of 10 days, or less depending on humidity level and temperature. Thereafter it must be disposed of. ● The treatment can be stopped and restarted any time by pressing the On/Off Button [7].
- **Taking off the Pressure Chamber [1]:** (a) Turn off the Control Unit [4] by pressing the On/Off Button [7]. (b) Roll down the Seal [20] around the neck of the Pressure Chamber [1] shell. (c) Push the Release Valve Button [28] to deflate the Padding [17]. (d) Gently remove the leg from the Pressure Chamber [1] after deflating the Padding [17]. (e) Roll up the Seal [20] from the Pressure Chamber neck to release it from tension. Do not leave the Seal [20] stretched over the Pressure Chamber neck.
- **After daily treatment:** If possible, leave the Pressure Chamber [1] and Control Unit [4] in place for daily treatment procedures. If the Control Unit [4] and Pressure Chamber [1] are left in place, make sure they are close to each other, to avoid tripping hazard, and unplugged from the electrical socket. ● If you used an Absorgel Pouch [26], seal it inside a sealable plastic bag until next use. ● If you need to move the system, use the Control Unit Handle [30] and the Pressure Chamber Handle [31] ● Wound dressings should be inspected after each treatment and replaced if necessary. ● If the Hose with Filter [24] is clogged, please contact your clinician or your local distributor.

- **Providing data to your clinician:** (a) Ensure that the Control Unit [4] is turned off but plugged into an electric socket (b) Insert the blue USB Memory Stick [32] supplied with the Pressure Chamber [1] into the USB Memory Stick Port [33] (c) Press the On/Off Button [7]. (d) When the Control Unit [4] has detected the USB Memory Stick [32], it will start to transfer saved treatment data onto it. The progress bar fills up from left to right and the Data Transfer symbol [34] appears on the Display [10]. (e) When the data transfer is complete, the Data Transfer symbol [34] and the Smiley Face symbol [9] are displayed. (f) Turn off the Control Unit [4] by pressing the On/Off Button [7] and then Remove the USB Memory Stick [32].
- **Adding time to system:** When the available time goes down to 25 hours or less, the Control Unit [4] will display an hourglass and the available time left on the Multi Screen [8]. Consult your clinician if you need a new Timecard [5]. When the time is up, an hourglass, the text "Oh" and a General Error symbol [35] will be displayed on the Display [10]. You may continue to use the Control Unit [4] when you have inserted a new Timecard [5] into the Control Unit [4]. To do so, follow the steps as described in section "Activate FlowOx™ with the Timecard."

VIEWING TREATMENT DATA (CLINICIAN)

Note: Connecting a PC to the Control Unit may only be done by a trained professional. PC must only be connected to the Control Unit outside of the patient area. Before connecting the Control Unit to a PC with the USB, make sure the PC has the latest security updates and an up-to-date antivirus software installed.

Installing FlowOx™ Software: FlowOx™ Software enables the clinician to view the treatment data, get or delete treatment data from the Control Unit [4], synchronize the Control Unit [4] time with the local time on the PC, and view the remaining time available on the Control Unit [4]. To install the software, (a) Insert the red USB stick with the FlowOx™ Software [36] supplied with the Control Unit [4] into the PC. (b) Copy the file "FlowOx-PC-SW" from the USB stick onto the PC. (c) Test the software by opening the file. **Note:** FlowOx™ will still work in cases where software cannot be used.

Installing the driver: The communication with the Control unit requires that the specific driver is installed on the PC where you are running the FlowOx 2 software. **Note:** You can read treatment data from the USB stick without installing the driver. Only the connection to the Control unit requires the driver installation. **Note:** To install the driver you shall have the administrator privileges. (a) Insert the red USB stick with FlowOx™ Software [36] into the PC. (b) The USB stick contains the file with the name: **InstallDriver**. Double click the **InstallDriver** file and follow the installation instructions (c) Restart your PC after installation. **Stored treatment data** can be viewed by (a) Launching the FlowOx™ Software (b) Clicking on the Folder icon [37] and then on the desired file. The file can be found on the USB Memory Stick [32] (see section "Providing data to your clinician") or in a previously created PC folder.

Direct Control Unit treatment data can be viewed in the following way: (a) Switch off the Control Unit [4]. (b) Connect the PC with the USB Cable [11] to the USB Cable Port [13] on the Control Unit [4]. (c) Turn on the Control Unit [4]. (d) Launch the FlowOx™ Software. (e) Click on the Cable icon [14] and then click on the Get Data icon [38]. • The treatment data that is stored on the Control Unit [4] can be deleted by pressing the Trashcan icon [39]. • The time left on Timecard is shown on the bottom right side "Control Unit Info" screen. • **Note:** The USB Cable [11] can be used in locations where portable memory devices are not permitted.

The treatment data is shown from 3 viewpoints (tabs), (a) The treatment summary shows the daily average treatment time, average daily session length, days used, and total hours used for a specific period. (b) The treatment days shows the accumulated daily usage divided into sessions for a specific period. (c) The treatment details show individual sessions per day, at the specific time they occurred, along with any errors. • In each of the viewpoints, it is possible to save data from a specified time period, save screen image or print. • To close the software, select "X" on the right upper corner.

SAFETY AND WARNINGS

FlowOx™ is intended for indoor use only. Do not use in rooms where the device could get wet such as in a bathroom. FlowOx™ should only be used in the conditions specified in Environmental and handling conditions • The FlowOx™ device should be kept away from direct sources of heat such as an open fire • Care should be taken to avoid tripping over the Hose with Filter [24] and the External Power Supply [2] • Caution should be taken to avoid accidental strangulation of yourself or others with inappropriate use of the Hose with Filter [24] and the power cord • Care should be taken with small parts to avoid accidental ingestion as this could cause suffocation • The power cable must never be attached to the wall power socket if it is disconnected from the Control Unit • Do not use with equipment not specified in the instruction manual, doing so could harm the operator • Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally • Portable RF communications equipment, such as mobile phones or other mobile devices

equipped with WIFI or Bluetooth, should not be used closer than 30 cm (12 inches) to any part of the Control Unit [4] or External Power Supply [2] and its cables • Use the system only with cables and other electrical accessories provided by the manufacturer. Using other cables or accessories could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation • Keep the FlowOx™ system away from children and pets • Make sure that the positioner is installed correctly inside the Pressure Chamber to avoid any pressure points on the foot • Make sure that there are no pressure points on the foot or leg from the Positioner, Seal, Padding, or Pressure Chamber • Do not try to run the Control Unit for more than 12 hours a day • The clinician and patient should be trained before using FlowOx™

- The Pressure Chamber [1] is a single patient use device. Do not share the Pressure Chamber [1] with different patients since this could lead to cross contamination • Do not rest or push the non-treated leg on the Pressure Chamber [1] shell during treatment. Doing so could cause damage to the skin surface of the area of the leg resting or pushing on the Pressure Chamber [1] • Be sure that your foot and leg are placed properly in the Pressure Chamber [1] without touching the Pressure Chamber [1] shell. Touching the inside surface could cause a pressure point, resulting in tissue injury or disturbance to wound dressings
- The inside of the Pressure Chamber [1] needs to be visually inspected both prior to and after use for any contamination such as wound exudate, dirt, and fibres • Clean the Pressure Chamber [1], Padding [17], Seal [20], and Positioner [16] in accordance with section “Cleaning”, if visible contamination is present • Check that the Hose with Filter [24] is visibly clean. Do not use the Pressure Chamber [1], if the Hose with Filter [24] contains any blood or liquid, in this case discard the Pressure Chamber [1] and use a new one • Care should be taken when putting on and taking off the Pressure Chamber [1], so as not to disturb wound dressings • Do not use with uncovered wounds • A non- shedding sock/tubular liner should be worn over all dressings and between the leg and the Positioner [16] and Padding [17] contact areas • Shedding socks should not be used as they potentially clog the pneumatic circuits and reduce in-use life • Do not attempt to modify the FlowOx™ Control Unit. The Control Unit cannot be serviced. When broken or malfunctioning, contact your clinician. Modifying the unit could harm the operator • Remain seated. Standing up or walking could lead to injury • Use of FlowOx™ has been reported in patient studies to be relaxing and injury may result if you are not adequately supported should you fall asleep • Do not attempt to lift FlowOx™ if you have reduced strength. It could fall out of your hands and injure you or someone else • Patients with a rapidly deteriorating condition should not start using FlowOx™ • FlowOx™ should be used in addition to standard of care in all cases • System parts must be verified before usage • In case of system failure, contact the clinician • USB should not be removed until device is correctly powered off • Correct disposal of the Control Unit is to dispose of it as electrical waste • Used Pressure Chambers should either be cleaned with hypochlorite solution prior to disposal as non-hazardous domestic waste or be disposed of as hazardous waste in a hospital or professional healthcare setting as per the National Waste Guidelines. Not doing so could risk spreading diseases to others • Visual inspection of the Absorgel Pouch [26] should be done prior to use to ensure that there is no rupture on it

CARE AND MAINTENANCE

Maintenance: The operator is not required to perform any preventive maintenance on the equipment. In the event of a fault, please contact your clinician or your distributor.

Cleaning: The FlowOx™ device can be cleaned if needed with wipes moistened with water and/or alcohol impregnated wipes.
Note: Do not get the Control Unit [4] wet • Do not use any abrasives or oxidizing agents such as detergents or hypochlorite (bleach) if the FlowOx™ system is to be used for further treatment.

Recyclability: Used Pressure Chambers should be cleaned with hypochlorite solution prior to disposal as per the National Waste Guidelines. Not doing so could risk spreading diseases to others • Do not dispose of the Control Unit [4] with your general waste. It contains a battery and other electronic components and must be disposed as electronic waste. When its end-of-life is reached, contact Otivio or your authorized distributor • The Pressure Chamber [1] is a disposable item and has a lifetime of one year.

Replacement parts: Please contact your clinician or distributor if you need replacement parts., i.e., Padding (900-00024) [17], Seal (900-00026) [20], 2x Absorgel Pouch (900-00025) [26], Refurbishment Kit (900-00028) containing 1x Padding [17] and 2x Seal [20] and/or the Strap (900-00029) [23]. These spare parts/consumables are not covered by the MDR 2017/745 certification.

TROUBLE SHOOTING

If an error occurs, the Treatment Ongoing symbol [27] disappears from the Display [10], Error code is shown on the Multi Screen [8] and a General Error symbol [35] or Air Leakage symbol [40] is shown. If an error code appears on the Multi Screen [8], it should be documented, and the clinician or distributor informed as soon as possible. Error Codes:

E1 Minor Air Leakage | Higher leakage than normal but treatment still being delivered effectively ● Solution: The clinician or distributor should advise patient at the next check-up, that it is important to adjust the Seal [20] and the Hose with Filter [24] so that there is minimal leakage. **Note:** this error is only shown in the data log

E2 Major Air Leakage | Prolonged leakage, treatment is not being delivered ● Solution: Check that the Seal [20] is sealed air tight to the bare leg, check that the Hose with Filter [24] is tightly connected to the hose connector [25], check that the air vents [19] are correctly inserted into the holes in the PCH and adjust if necessary. Restart the device by pressing On/Off Button [7]

E3 Low Air Volume | Not enough air to function (Low air volume) ● Solution: check for obstruction, e.g., kinked hose. Restart the device by pressing On/Off Button [7]

E4 Data transfer error | Data not transferred, e.g., due to taking out USB Memory Stick [32] too early or software error ● Solution: Switch off the device and try inserting the USB Memory Stick [32] into the Control Unit [4] again ● Note: if the problem persists, contact the clinician or local distributor ● The Control Unit [4] can be reset by pressing On/Off button [7] for 13 seconds.

COMPLAINTS AND ADVERSE EVENTS

Complaints about FlowOx™, including any adverse events experienced in its use, should be communicated to your distributor or directly to Otivio at info@otivio.com. Serious adverse events should also be reported directly to Otivio at info@otivio.com and the relevant authorities where you live, contacts are shown below:

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Norway: meddev-no@legemiddelverket.no

SPECIFICATIONS

Classification of the FlowOx™ system: Class IIa in accordance with the Medical Devices Regulation (MDR 2017/745) ● Class A in accordance with IEC 62304:2006+A1:2015, Software Safety Classification

Implemented Safety Standards: IEC 60601-1:2005+AMD1:2012+AMD2:2020 Ed. 3.2 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance ● IEC 60601-1-2:2014+AMD1:2020 Ed.4.1 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance. Collateral Standard: Electromagnetic disturbances - Requirements and tests ● IEC 60601-1-6:2010+AMD1:2013+AMD2:2020 Ed. 3.2 Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance. Collateral standard: Usability ● IEC 60601-1-11:2015+AMD1:2020 Ed. 2.1 General requirements for basic safety and essential performance. Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment ● IEC 62304:2006+A1:2015 Medical device software – Software life-cycle processes ● BS EN ISO 14971:2019+A11:2021 Medical devices – Application of risk management to medical devices ● BS EN ISO 10993:2020 Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process ● IEC 62366-1:2015+AMD1:2020 Ed. 1.1 Medical devices – Application of usability engineering to medical devices ● BS EN ISO 15223-1: 2021 Medical devices - Symbols to be used with

medical device labels, labelling and information to be supplied. Part 1: General requirements • ISO 15223-2 Rev 1.0: 2010 Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied. Part 2 • BS EN ISO 20417: 2021 Medical devices — Information to be supplied by the manufacturer • EN 15986:2011 Symbol for use in the labelling of medical devices - Requirements for labelling of medical devices containing phthalates

Manufacturing Standards used: EN ISO 13485:2016+A11:2021 Medical devices – Quality management systems – Requirements for regulatory purposes.

Environmental and handling conditions: Operation Temperature range: +15 to +40°C. **Note:** if device is taken from storage that is outside of operating temperature (15 – 40°C), wait at least 2 hours before using the device.

• Operation Air Humidity: 15 – 95 % • Operation Ambient pressure: 70 – 106 kPa (525 – 795 mmHg) • Storage and transport Temperature range: - 25 to +70 °C • Storage and transport Air Humidity: Up to 93 % relative humidity at +70 °C • Storage and transport Ambient pressure: 50 – 106 kPa • Degree of enclosure protection: IP21 (Protected against solid objects over 12mm e.g., fingers. Protection against vertical liquid drops) • Vibration/Shock/Bump: It is possible to transport the system world- wide by air, road, ship and train • Drop/Free fall: It is possible to transport the system world- wide by air, road, ship and train • EMC/ESD: The FlowOx™ system meets the requirements in accordance with IEC 60601-1-2 Electromagnetic compatibility • Life expectancy: Control Unit [4]: 3 years, Pressure Chamber [1]: 1 year, Padding [17]: 3 months, Seal [20]: 3 months, Strap [23]: 1 year.

Note: The pressure cycle of the FlowOx™ system as specified in section “How It Works” may be affected if the system is used in an environment with excessive electromagnetic interference. Excessive interference may cause erratic behavior and/or inaccurate vacuum and timing of the treatment cycle. Please see section “EMC Specification” for guidance on electromagnetic interference and appropriate conditions for use.

Technical specifications: Description: Enhancement of blood flow to lower limbs • Parts/Names: FlowOx™ 2.0 (900-00012), FlowOx™ 2.0 Control Unit (900-00017), FlowOx™ 2.0 Firmware version 1.0.9 (510-00006), FlowOx™ 2.0 Pressure Chamber (900-00018), FlowOx™ 2.0 PC SW version 1.0.7 (520-00003), Mains power cord, Instruction manuals • Control Unit dimensions: 340 x 250 x 200 mm • Pressure Chamber size: Max EU 46 / UK 11½ • The Pressure Chamber will fit all patients with leg circumference of minimum of 25cm and maximum of 47cm, measured at 20cm above the lateral malleolus • Weight: Control Unit: 3 kg, Pressure Chamber: 3 kg • Pressure accuracy: ±5mmHg • Operation time: Normal 2h/day • Power supply: 100-240 V~, 50-60 Hz, 1,4 – 0,7 A • Power consumption: Max 60 VA • Battery: Lithium cell battery, CR2032, 230MAH, 3V • Computer requirements: Min. 1.5 Giga-byte RAM Windows 7 or later • Min. resolution: 800x600pix. Min. 24-bit colour support

DEFINITION OF SYMBOLS ON THE DEVICE AND PACKAGING

Refer to the symbols from the cover page of this instruction manual.

Display: [9] Smiley Face symbol (indicates that an action has been completed) • [12] Clock Set symbol (indicates that the Control Unit [4] time is being set) • [27] Treatment Ongoing symbol (indicates that the treatment is ongoing) • [34] Data Transfer symbol (indicates the transfer of data from Control Unit [4] to the USB Memory Stick [32]) • [35] General Error symbol (indicates that an error has occurred) • [40] Air Leakage symbol (indicates that air is leaking somewhere in the system)

PC software: [14] Cable icon • [15] Clock Set icon • [37] Folder icon • [38] Get Data icon • [39] Trashcan icon

On device and packaging: [41] Refer to the instruction manual for power supply connection. Refer to the instruction manual to ensure correct use of the device • [42] Keep dry • [43] Manufacturer • [44] Batch number • [45] Type reference • [46] Serial number, including manufacturing date • [47] Do not dispose as general waste • [48] Recyclable material • [49] Warning: Do not open, causes serious eye irritation (Absorgel Pouches [26]) • [50] Temperature transport conditions • [51] Type BF applied part • [52] Humidity transport conditions • [53] IEC Protection class II • [54] The Device Complies with Medical Device Regulation 2017/745 • [55] Warning: Dangerous voltage • [56] General warning sign • [57] Left Foot symbol (Indicates how to adjust the Positioner [16] for left foot usage. Positioner [16] should be placed in such a way so that the foot symbol is aligned to the left foot of the patient) • [58] Right Foot symbol (Indicates how to adjust the Positioner [16] for right foot usage. Positioner [16] should be placed in such a way so that the foot symbol is aligned to the right foot of the patient) • [59] Single patient multiple use • [60] Patient name • [61] Medical device • [62] Unique device identifier • [63] Date of manufacture • [64] Quantity • [65] Expiration Date

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EMC SPECIFICATION

Guidance and manufacturer's declaration – electromagnetic immunity		
The Otivio AS FlowOx™ unit is intended for use in the electromagnetic environment below. The customer or the end user of the FlowOx™ should assure that it is used in such an environment.		
Immunity test	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	
Electrical fast transient/burst IEC61000-4-4	± 2 kV for power supply lines 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0,5 kV, ± 1 kV Line-to-line	Mains power quality should be that of a typical commercial or hospital environment.

Voltage dips IEC 61000-4-11	0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the FlowOx™ requires continued operation during power mains interruptions, it is recommended that the FlowOx™ be powered from an uninterruptible power supply or a battery.
Voltage interruptions IEC 61000-4-11	0 % UT; 250/300 cycle (250 periods for 50 Hz or 300 periods for 60 Hz)	
Rated power frequency magnetic field IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	Minimum distance between the The Otivio AS FlowOx™ unit and sources of power frequency magnetic field of at least 15 cm
NOTE: Ut is the AC mains minimum and maximum rated input voltage		

Guidance and manufacturer's declaration – electromagnetic emissions		
The Otivio AS FlowOx™ unit is intended for use in the electromagnetic environment below. The customer or the end user of the FlowOx™ should assure that it is used in such an environment.		
Emissions tests	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group I	The FlowOx™ uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The FlowOx™ is suitable for use in all establishments, including domestic establishments.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity
The Otivio AS FlowOx™ unit is intended for use in the electromagnetic environment below. The customer or the end user of the FlowOx™ should assure that it is used in such an environment.
Electromagnetic environment - guidance
Portable and mobile RF communications equipment should be used no closer to any part of the FlowOx™, including cables, than the recommended separation distance calculated from the equation applicable to frequency of the transmitter.

Recommended separation distance

$$d = 0,6 * \sqrt{P}$$

where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).

Field strengths from fixed RF transmitters, as determined by an electro-magnetic site survey,^a should be less than the compliance level in each frequency range.

Interference may occur in the vicinity of equipment marked with the following symbol:

Immunity test	IEC 60601 Test level
---------------	----------------------

Conducted RF IEC 61000-4-6	3 Vrms 15 kHz to 80 MHz 6 Vrms 15 kHz to 80 MHz in ISM and amateur radio bands 80 % AM at 1 kHz
-------------------------------	---

Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2,7 GHz 80 % AM at 1 kHz
------------------------------	---

Proximity fields from RF wireless communications equipment IEC 61000-4-3	TEST FREQUENCY (MHz)	BAND ^{a)} (MHz)	SERVICE ^{a)}	MODULATION ^{b)}	MAXIMUM POWER (W)	DISTANCE (m)	IMMUNITY TEST LEVEL (V/m)
	385	380 - 390	TETRA 400	Pulse modulation ^{b)} 18 Hz	1,8	0,3	27
	450	430 - 470	GMRS 460, FRS 460	FM ^{c)} ± 5 kHz deviation 1 kHz sine	2	0,3	28
	710	704 - 787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9
	745						
	780						
	810	800 - 960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	2	0,3	28
	870						
	930						
	1720	1700 - 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^{b)} 217 Hz	2	0,3	28
	1845						
	1970						
	2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	2	0,3	28
	5240	5100 - 5600	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9
	5500						
5785							
NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.							
a) For some services, only the uplink frequencies							
b) The carrier shall be modulated using a 50 % duty cycle square wave signal							
c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.							

Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the FlowOx™ is used exceeds the applicable RF compliance level above, the FlowOx™ should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the FlowOx™.

Recommended separation distances between portable and mobile RF communications equipment and the FlowOx™

The FlowOx™ is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the FlowOx™ can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the FlowOx™ as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	380 MHz to 390 MHz $d=6/27*\sqrt{P}$	430 MHz to 470 MHz 800 MHz to 2,57 GHz $d=6/28*\sqrt{P}$	704 MHz to 787 MHz 5,1 GHz to 5,8 GHz $d=6/9*\sqrt{P}$
0,01	0,03	0,03	0,07
0,1	0,08	0,07	0,21
0,2	0,10	0,10	0,30
1	0,23	0,22	-
1,8	0,30	0,29	-
2	-	0,30	-

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where E is the immunity test level (V/m) and P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Manufacturer:

Otivio AS
Drammensveien 130
NO-0277 Oslo
Norway
www.flowox.com
info@otivio.com

CE 2460